Sedens 400

User Manual

Be well, Live well
Wellell

IMPORTANT SAFEGUARDS
READ CAREFULLY ALL INSTRUCTIONS BEFORE USING

Warning

- (1). Any serious incident that has occurred to a user or patient in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or patient is established.
- (2). Use this product only for its intended use as described in this manual. Do not use other cushions not recommended by the manufacturer.

NOTE, CAUTION AND WARNING STATEMENTS

NOTE: Indicates some tips.

⚠ **CAUTION:** Indicates correct operating or maintenance procedures to prevent damage to or destruction of the equipment or other property.
WARNING: Calls attention to a potential danger requiring correct procedures or practices to prevent personal injury.

SYMBOLS

	Refer to instruction manual/ booklet		Humidity Limitation
	Do Not Dry Clean		Do Not Iron
	Tumble Dry, Normal, Low Heat		Machine wash, regular/ normal, 95°C (203°F)
	Do Not Bleach		Temperature Limitation
	Manufacturer		Authorised representative in the European community
	Maximum/minimum user weight		Date of manufacture
	Keep Dry		Fragile, Handle with Care
	This Way Up		Use No Hand Hooks
	CE Marking		Medical Device
	Unique Device Identifier		Importer
	Constructed from recyclable materials		Do not use cutter

1. Intended Use

The device is intended to help prevent pressure injuries in patients at risk and to provide postural support for those who need stability. It can only be operated by personnel with adequate training in seating and positioning.

The device is suitable for use in acute care settings, long-term care facilities, and home environments.

1.1. Indication for use

The device is indicated for patients whose weights are within the specified minimum and maximum support weight limits.


1.2. Contraindication (when used with a pump)

Do not use on patients with unstable fractures, including unstable spinal fractures, or those undergoing cervical traction or skeletal traction.

2. Cushion Description

- (1) Stretch Cover
- (2) Zipper
- (3) Air cells with foam inside
- (4) Quick Connector (to Wellell pump)
- (5) Cap with Pressure Adjustment Valve


NOTE: Please refer to the pump manual for the respective pump description.

3. Installation & Adjustment


3.2. Alternating Hybrid Cushion Use- With Wellell Pump

NOTE: The Sedens 400 Hybrid Cushion is compatible with all pumps that use 1-in-2 Wellell Quick Connectors, except Auto systems.

- (1) Place the Sedens 400 hybrid cushion in the wheelchair seat so the user will be properly supported and immersed when in their normal sitting position.
- (2) Uncap the quick connector and connect to a compatible Wellell pump. A "click" sound will be heard when the connection is completed.
- (3) Install the pump according to the pump manual.
- (4) Have the user sit on the cushion. Check to make sure the cushion is properly sized to the user and wheelchair.
- (5) Adjust an appropriate pressure setting to ensure the patient's comfort without bottoming out.
- (6) Keep hand under user's lowest bony prominence and make sure there is between 1/2 in. (1.3cm) to 1 inch (2.5 cm) of cushion between the user's lowest bony prominence and the seating surface. If there is not enough space in between, please adjust the pump pressure to a higher level.

3.3. Temporary Static Cushion Use (without pump)

- (1) Capped on the quick connector and open air valve by turning the valve **counterclockwise**.
WARNING: Make sure the cushion is unweighted to inflate. (No patient should be sitting on the cushion at this time.)
- (2) When the cushion is fully inflated, close the air valve by turning the valve **clockwise**.
- (3) Have the user sit on the cushion. Check to make sure the cushion is properly sized to the user and wheelchair.
- (4) Turn the valve **counterclockwise** to release air, while keeping your hand under the user's lowest bony prominence. Release air until there is between 1/2 in. (1.3cm) to 1 inch (2.5 cm) of cushion between the user's lowest bony prominence and the seating surface.
- (5) Turn the valve **clockwise** to close.

4. Cleaning

Cover Material: Stretch	
-------------------------	--

Wipe down the cushion with a damp cloth pre-soaked with warm water containing a mild detergent, or chlorine bleach, followed by an approved intermediate level disinfectant. The top cover can also be completely removed for laundry with the indicated temperature; however, it is recommended that the user still check with local policy to determine the time/ temperature ratio required to achieve thermal disinfection. After cleaning, please avoid direct sunlight or dusty areas. All parts should be air-dried thoroughly before use.

- ⚠ **CAUTION:** Do not use phenolic-based products for cleaning.
- ⚠ **CAUTION:** Keep the cushion out of direct exposure to sunlight, and make sure to dry the cushion after cleaning.
- ⚠ **CAUTION:** Always perform the cleaning procedure before each new patient uses the cushion.

5. Technical Specifications

Cushion		Specification
Environment	Temperature	Operation: 50°F to 104°F (10°C to 40°C) Storage: 5°F to 122°F (-15°C to 50°C) Shipping: 5°F to 158°F (-15°C to 70°C)
	Humidity	Operation: 10% to 90% non-condensing Storage: 10% to 90% non-condensing Shipping: 10 % to 90% non-condensing
Dimensions (L x W x H)		44 x 41 x 9 cm / 17.3" x 16.1" x 3.5"
Weight		1.1 kg / 2.4 lb
Max. Weight Capacity		250 kg / 550 lb

NOTE:

- Pump specification/EMC information please refer to the pump manual.
- Consult the distributor for other technical documents.
- Please follow national requirements to dispose of the unit properly.
- The manufacturer reserves the right to change product specifications without prior notice.

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All rights reserved.
Wellell Inc.
Apex Medical (Kunshan) Corp.