Optima Turn

User Manual

Be well, Live well
Wellell

CONTENT

ENGLISH

USER MANUAL 1

ESPAÑOL

MANUAL DE INSTRUCCIONES 20

PORTUGUÊS

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FRANÇAIS

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DEUTSCH

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DUTCH

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ITALIANO

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MODEL NO.: 9P-085500

PLEASE READ THE MANUAL BEFORE USE

IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE USING

DANGER - To reduce the risk of electrocution:

- (1). Always unplug this product immediately after using.
- (2). Do not use while bathing.
- (3). Do not place or store this product where it can fall or be pulled into a tub or sink.
- (4). Do not place in or drop into water or other liquid.
- (5). Do not reach for a product that has fallen into water. Unplug immediately.

WARNING- To reduce the risk of burns, electrocution, fire, or injury to persons:

- (1). Evaluate patients for entrapment risk according to protocol and monitor patients appropriately.
- (2). The product may be used with caution for patients with spinal injury, but suggested to consult with physician before use. However, it should not be used for patients with unstable spinal fractures.
- (3). Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may result from a child swallowing a small part detached from the device.
- (4). Use this product only for its intended use as described in this manual. Do not use other mattress not recommended by the manufacturer.
- (5). Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to your supplier or Wellell for examination and repair.
- (6). Keep the cord away from heated surfaces.
- (7). Never block any air openings of this product or place it on soft surfaces, such as a bed or couch, where openings may be blocked. Keep the air opening free of lint, hair, and other similar particles.
- (8). Never drop or insert any object into any opening or hose.
- (9). Do not modify this equipment without authorization of the manufacturer.
- (10). Mattress covers have passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.
- (11). Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.
- (12). Third conductor in the POWER SUPPLY CORD is only a functional earth.
- (13). The Operator cannot operate the instrument and touch the patient at the same time.
- (14). If the RTC battery needs to be replaced, please contact local representative or place of purchase.

CAUTION-

If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.


NOTE, CAUTION AND WARNING STATEMENTS:

NOTE- Indicate some tips.

CAUTION- Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

WARNING- Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

SYMBOLS

	<p>Authorized representative in the European community.</p>		<p>Temperature Limitation</p>
	<p>Manufacturer</p>		<p>Dry clean, Any Solvent Except Trichloroethylene</p>
	<p>Complies with standards protecting against electric shock for type BF equipment.</p>		<p>Do Not Iron</p>
	<p>Consult operating instructions for use</p>		<p>Tumble Dry, Normal, Low Heat</p>
	<p>Attention, should read the instructions.</p>		<p>Do Not Tumble Dry</p>
	<p>Attention, should read the instructions.</p>		<p>Do Not Bleach</p>
	<p>Functional earth. (for UL only)</p>		<p>Do Not Dry Clean</p>
	<p>Class II</p>		<p>Machine wash, regular / normal, 95 degrees C (203 degrees F)</p>
<p>IP21 </p>	<p>Protected against solid foreign objects of 12.5 mm and greater; Protection against vertically falling water drops</p>		
	<p>Attention – Observe proper Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.</p>		

1. Introduction

This manual should be used for initial set up of the system and for reference purposes.

1.1. General Information

The system is a high quality mattress system suitable for treatment and prevention of pressure injury. The system has been tested and successfully approved to the following standards:

IEC/EN 60601-1
IEC/EN 60601-1-2
IEC/EN 61000-3-2 Class A
IEC/EN 61000-3-3
CISPR 11 Group 1, Class B

EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

1.2. Intended Use

This product is intended to help and reduce the incidence of pressure injury while optimizing patient comfort. It also provides following purposes:

- for long term home care of patients suffering from pressure injury.
- for pain management as prescribed by a physician.

The product can only be operated by personnel who are qualified to perform general nursing procedures and have received adequate training in knowledge of prevention and treatment of pressure injury.


NOTE- Equipment is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.

2. Product Description

Unpack the box to check if there are any damage which may have occurred during shipment. Please contact local representative or place of purchase.

2.1. Pump and Mattress System


- ① CPR
- ② Pump unit
- ③ Mattress system
- ④ Quick connector

2.2. Pump Unit


Front

- ① Quick connector
- ② Front panel


Rear

- ① Power switch
- ② Air filter
- ③ Fuses
- ④ Mounting brackets

2.3. Front Panel

(1). Power Key 

Power Key is to turn on and turn off pump.

(2). Pause Key 

- Allow caregivers to temporarily stop the action when emergency occurred. Caregivers could change to other settings during this period of time which will return to the previous mode automatically after 10 min.
- To go to the previous mode setting please press the pause key again.

(3). Panel Lock Key 

User interface will be locked automatically 5 minutes after last operation or by pressing Panel Lock key. To unlock interface, press and hold Panel Lock key for 1 second.

NOTE- Always unlock user interface before turning off the pump.

(4). Max Firm Key 

The surface will go into Max firm automatically when the power is started.

NOTE- The LED light of Max Firm mode will be on constantly until the mattress is reached to the maximum pressure setting.

NOTE- This key is to make sure mattress can be reached maximal operating pressure. This mode will last for 20 minutes and then return to previous setting automatically. This function can be used during patient transfers or nursing procedures for better support. To disable this function, please press Max firm key again.

(5). Modes A. Alternate Mode 

The air cells inside the mattress are deflated and inflated periodically to avoid prolonged pressure in certain areas of patient.

B. Continuous Low Pressure 

Function LED indicator is steady-on in constantly non-alternating mode. The system will redistribute body mass over a greater surface area at a constant low pressure.

C. Seat Inflation 

Press the  key to enter the desired mode in alternate or CLP mode for additional supports to the patient during upright position without bottoming out.

(6). Turning Mode 

Three turning modes including  LEFT TURN,  RIGHT TURN and  LEFT-SUPINE-RIGHT TURN can be selected. Turning angle is from patients' direction and hold at that position for dwell time. Top air cells are non-alternating during the turning process, and will change to alternating mode when it reaches the position of the desired angle and when it comes back to supine position as well.

NOTE- When dwell time of 10 min is selected, even reaches the position of the desired angle it doesn't change to alternating mode.

(7). Turning Angle  (Up to 30 degree turning under Full-turn Mode, subject 75Kg with shoulder width 45.7cm)

-  Light-turn,  Half-turn,  Full-turn
- For One-Time Turning which indicated with LED flashing, press and hold  key for 3 seconds. When turning is completed, it will go back to supine position with alternating mode.

(8). Pressure Adjust Levels

- Press the / to get the correct setting of the pressure level based on patient weight.
- The mattress pressure can be manually adjusted for patient comfort using the key  for a firmer setting and decrease  for a softer setting.

NOTE- When the power is started (i.e., the surface is inflated), it will automatically go to "Max Firm" mode to hasten inflation. Once the system is ready, the system will enter CLP mode for 10 minutes. If there are no changes of modes, system will switch to Alternate Mode.

NOTE- You can also check if the pressure is suitable for the patient by sliding one hand beneath the air cells at the level of the patient's buttocks. Always leave at least 1" inch space between patient and the static cell to prevent bottoming out.

(9). Cycle Time 

- By pressing this button after changing the operation mode, **Alternate Mode Cycle Time** can be selected between 10, 15, 20 and 25 minutes.
- By pressing this button after changing to the turning mode, Turning **Dwell Time** can be selected between 10, 20, 30 and 60 minutes.

NOTE- Cycle time and dwell time duration can be selected based upon patient comfort and desired outcome.

(10). Alarm Mute Key 

When there is an alarmed situation happened, the number display will show corresponding error code with buzzer sounded every few seconds. The buzzer can be temporarily suspended by pressing the **Alarm Mute Key**. Alarm buzzer will signal repeatedly every 3 minutes until the alarmed issue is solved.

A. PFA (Power Failure Alarm) 

The LED light and buzzer will activate when the power failure happens. To disable both buzzer and LED, press the mute button.

NOTE- The pump can seal the air inside the mattress, which temporarily keeps the patients from bottoming out during power failure. Patients should be moved if the power failure issue cannot be solved soon.

B. Low Pressure Indicator 

When low pressure LED lights up, the pressure inside the air mattress is below normal range. Please refer to troubleshooting guide.

NOTE- The pump can seal the air inside the mattress, which temporarily keeps the patients from bottoming out during low pressure condition. Patient should be moved if the power failure issue cannot be solved soon.

C. Tech. Support 


The Service indicator light will activate when an issue occurs related to the system. Please contact the local representative or place of purchase for technical support.

3. Installation

3.1. Pump & Mattress Installation

 NOTE- Unpack the box. Please contact your support personnel or package vendor if there are any damage occurred during shipment.

 NOTE- USB function is provided for use by distributors or authorized technicians.

(1). Place the mattress on top of bed frame and foot print toward foot end. Secure the mattress firmly by fixing the straps to bed frame to ensure mattress will be free from move. Securing straps is located on base of the mattress

 NOTE- Make sure that ALL the Heel Relief Connectors and emergency CPR knob are closed and locked in position.


 NOTE- Make sure the bed side rails are raised.

(2). Hang the pump onto bed rail (foot-end), and adjust hangers to best upright position of the pump, or place the pump on flat surface.

(3). Connect air hose connectors of air mattress to pump unit. When a "click" sound is felt or heard, the connection is completed and secured.

 NOTE- Check and ensure the air hoses are not kinked or tucked under mattress.

(4). Plug the power cord into electrical outlet.

 NOTE- Make sure the pump unit is suitable for the local power voltage.

 NOTE- The plug is also served to disconnect the device.

(5). Then turn on power switch to ON position.

 CAUTION- The pump can only be applied to the mattress recommended by the manufacturer. Do not use it for any other purpose. (Applied part: air mattress) Do not position the equipment so that it is difficult to operate the disconnecting device.

- (6). The cable management tape is provided along the side of the mattress to manage power cable. Please make sure the wire has no signs of damage. After installation, extra length of the power cord, if any, should be neatly managed to prevent stumble accident and should be clear of moving bed mechanisms or other possible entrapment area.

NOTE- The EQUIPMENT should be firmly placed at position where caregivers can access easily.

4. Operation

NOTE- Always read the operating instruction before use.

4.1. General Operation

- (1). Switch the power on, the mattress starts the initial inflation. The Maxfirm light flashes until the mattress is fully inflated which takes approximately 30 min.
- (2). When the inflation is completed, the system will enter CLP mode for 10 minutes. If there are no changes of modes, system will switch to Alternate Mode.
- (3). Use the Pressure Adjustment key to get the correct setting of therapeutic pressure level based on patient weight.

NOTE: Pressure levels with the weight guide on the panel is only an indication and clinical judgement should always be taken into consideration.

- (4). For lateral turning, choose the Left-SUPINE-RIGHT turn key to get standard turning which helps the patient to turn to both sides and back to supine position, or either one direction in Left turn or Right turn mode.


NOTE: Do the hand check to check out if the pressure reduction is adequate and the patient has not bottomed out.

NOTE: For patient comfort, always use the head pillow when patient is under turning mode.

4.2. Emergency CPR Operation

NOTE- When there is a need to perform CPR on the patient, quickly turn the CPR knob to release air from mattress. CPR knob is located at the right-head side of the mattress. The quick connector on the pump unit can be disconnected for even faster deflation.

4.3. Heel Relief Function

The last five cells from the foot-end have heel-relief connectors. Carers can easily disconnect the connector and heels can be suspended in air to completely eliminate pressure exerted on the heels.

For permanent off-loading/pressure relief:

- Select **no more than two** cells directly under the area that desired to relieve.


NOTE- This one to two cells can be left permanently deflated.

NOTE- Deflating more than two adjacent cells may affect the support of the patient during the normal alternating cycle.

4.4. Transportation

NOTE- When there is a need of patient transportation, pump with the power cord can be disconnected from the mattress during transport. For optimum comfort, change the system to CLP Mode and wait for 5 minutes before disconnecting is recommended. Cap on the quick connector as shown on the left to keep the air inside the mattress. The mattress remains inflated at least 24 hours depending on the setting.

5. Cleaning

It is important to follow the cleaning procedures to avoid cross contamination. Be sure to clean the surface in a dry and dust free environment. Wipe down the pump unit with a damp cloth pre-soaked with a mild detergent. Avoid contact with dust and proximity to dusty areas. Make sure that any cleaning agents you use will not harm or corrode the plastic casing on the pump unit. If your doctor or medical facilities have other special cleaning instruction, please follow the professional instruction.

CAUTION- Do not immerse or soak pump unit.

WARNING- Do not remove the casing of the pump to avoid electrical shock. All disassembly or repair should be done by authorized technicians.

CAUTION- The pump does not need oil lubrication; please do not disassemble the system.

<p>Cover Material: Stretch Silver* Stretch</p>	
--	--

Wipe down the mattress with 70 % alcohol solution or chlorine solution concentration must be less than or equal to 1,000 ppm. After disinfecting the mattress, use water and a soft cloth to wipe off all traces of disinfectant and make sure the mattress is dry before storage. Also the mattress top cover can be completely removed for laundry with water temperature indicated on the laundry label; however, it is recommended that user should check with local policy or regulation to determine the time/temperature ratio required achieving thermal disinfection. The cover may also be cleaned using sodium hypochlorite diluted in water. After cleaning, please avoid dust and proximity to dusty areas and all parts should be air dried thoroughly before use.

CAUTION- Do not use phenolic based products for cleaning.

CAUTION- After cleaning, dry the mattress without direct exposure of sunlight.

CAUTION- Do not use ozone or UV disinfection as these methods can damage the surface coating.

6. Storage

- (1). To store the mattress, lay the mattress out flat and upsides down.
- (2). Roll from the head end towards the foot end with CPR valve open.
- (3). Foot-end strap can then be stretched around the rolled mattress to prevent unrolling.

NOTE- Do not kink, crease or stack the mattresses and do not store the system in direct sunlight, high temperature or moisture area.


7. Maintenance

7.1. General

- (1). Check main power cord and plug if there is abrasions or excessive wear.
- (2). Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are stubbed together correctly.
- (3). Check airflow from the opening of the tubes in the quick connector. The airflow should alternately blow every half of the cycle time between two tubes at the lower right corner and one of the the other two tubes at the lower left corner will blow to give desired turning angle when it is under turning mode.
- (4). Check the air hoses if there is kink or break. For replacement, please contact local representative or place of purchase.

7.2. Fuse Replacement

- (1). Disconnect the plug from mains power when a blown fuse is suspected.
- (2). Remove the cover of the fuse holder by means of a flat screwdriver.
- (3). Insert a new fuse of the correct rating in, and replace the cover of the fuse holder back. The fuse should be rated as T1AL/250V type and VDE approved.

7.3. Air Filter Replacement

- (1). Open the air filter cap located at the back of pump.
- (2). Replace with a clean filter and put the cover lid back. The filter is reusable and can be washed gently with a mild detergent and water. Dry the filter before use.
- (3). Suggested to check and replace air filter regularly every six month.

8. Expected Service Life

The products are intended to offer safe and reliable operation when use or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indication on products. Otherwise, service and inspection of the devices generally should not be required.

9. Trouble Shooting

Q1 Power is not ON

- Check if the plug is connected to mains.
- Check for a blown fuse.

Q2 Low Pressure Alarm is on

- Check if the Quick Connector is tightly secured.
- Check if all tubing connections along mattress are secured.
- Check if the air hoses are kinked or broken.
- Check if the CPR knob is tightly closed.

Q3 Power Failure Alarm is on

- Check if the power is suddenly shut down.
- Check if the power cord is connected properly.

Q4 Patient is bottoming out

- Pressure setting might be inadequate for the patient, increase comfort range 1 level higher and wait for a few minutes for best comfort.

Q5 Mattress form is loose

- Check if all the snap buttons or straps of mattress are all securely fastened.
- Check if the mattress is fixed to the bed frame by straps.

Q6 No air released from some air outlets of the air tube connector

- This is in alternating mode. The air is released through outlets by turns at each cycle.

If the problems you have are not solved by the above information, please contact local representative or place of purchase.

10. Technical Specifications

Item		Specification
Pump		
Power Supply (Note: See rating label on the product)		AC 230V 50 Hz, 0.25A (for 230V system)
Fuse Rating		T1AL, 250V
Cycle time (Alternating)		10, 15, 20, 25 min.
Dwell time (Turning)		10, 20, 30, 60 min.
Dimension (L x W x H)		34.1 x 16.5 x 26.0 cm
Weight		4.5 Kg / 9.9 lb
Environment	Atmospheric Pressure	700 hPa to 1060 hPa
	Temperature	Operation: 10°C to 40°C (50°F to 104°F) Storage: -15°C to 50°C (5°F to 122°F) Shipping: -15°C to 70°C (5°F to 158°F)
	Humidity	Operation: 10% to 90% non-condensing Storage: 10% to 90% non-condensing Shipping: 10% to 90% non-condensing
Classification		Class II, Type BF, IP21 Applied Part: Air Mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)
Mattress		
Model		5" cells & 10" side bolster
Dimension (L x W x H)		200x85/90x13cm 78.7"x33.5/35.4"x5"
Weight		12.5 Kg / 27.5 lb
Max. Support Weight		180 Kg / 397 lb

NOTE-

- Consult the distributor or EU representative for further technical documents.
- The specification is also suitable for other areas operating with same power supply.
- The manufacturer reserves the right to modify the specifications without notice.

11. Error Code Table

Code No.	Description
00	Initial inflation overtime
01-06	Low pressure
08-12	System malfunction
16	Pressure control system malfunction
17	System access control malfunction
18	Mattress disconnected
19	Connect to the wrong mattress

Appendix A: EMC Information

Guidance and Manufacturer’s Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	
<p> Warning:</p> <ol style="list-style-type: none"> 1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used. 2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. 3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. 		

Guidance and Manufacturer’s Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	±8 kV contact ± 2, 4, 8, 15 kV air		±8 kV contact ± 2, 4, 8, 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2 kV, 100 kHz for power supply line ±1 kV, 100 kHz for input/output line		±2 kV, 100 kHz for power supply line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 0.5, 1 kV (line to line) ± 0.5, 1, 2 kV (line to earth)		± 0.5, 1 kV (line to line)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		230 Vac	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from

			80 MHz 80 % AM at 1 kHz	the equation applicable to the frequency of the transmitter.
Radiated RF EM Fields IEC61000-4-3	3 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz 380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation, (1kHz sine)	10 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz 380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation, (1kHz sine)	10 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz 380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation, (1kHz sine)	Recommended separation distance $d = \sqrt{P}$ 150kHz to 80MHz $d = 0.6\sqrt{P}$ 80MHz to 800MHz $d = 1.2\sqrt{P}$ 800 MHz to 2.7GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: U _T is the a.c. mains voltage prior to the application of the test level NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people				

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \sqrt{P}$	80 MHz to 800 MHz $d = 0.6\sqrt{P}$	800 MHz to 2.7 GHz $d = 1.2\sqrt{P}$
0.01	0.1	0.06	0.12
0.1	0.31	0.19	0.38
1	1	0.6	1.2
10	3.1	1.9	3.8
100	10	6	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Print-2022/All rights reserved
Apex Medical Corp.