ACS TURN

User's Manual

Care for a Healthy Life

CONTENTS

ENGLISH

USER MANUAL 1

ESPAÑOL

MANUAL DE INSTRUCCIONES 20

PORTUGUÊS

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FRANÇAIS

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DEUTSCH

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ITALIANO

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DUTCH

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Model No.: 9P-060120

PLEASE READ ALL INSTRUCTIONS BEFORE USE.

IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE OPERATING THIS DEVICE.

DANGER - To reduce the risk of electrocution:

1. Always unplug this product immediately after using.
2. Do not use while bathing.
3. Do not place or store product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquid.
5. Do not reach for a product that has fallen into water. Unplug immediately.

WARNING - To reduce the risk of burns, electrocution, fire, or injury to persons:

1. This product should never be left unattended when plugged in.
2. Evaluate patients for entrapment risk according to facility protocol and monitor patients appropriately. Choking accident may result from a child swallowing a small part that has become detached from the device or its accessories.
3. Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may result from a child swallowing a small part detached from the device.
4. Use this product only for its intended use as described in this manual. Do not use other mattress not recommended by the manufacturer.
5. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to a service center for examination and repair.
6. Keep the cord away from heated surfaces.
7. Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where their openings may be blocked. Keep the air opening free of lint, hair, and other similar particles.
8. Never drop or insert any object into any opening or hose.
9. Connect this product to a properly grounded outlet only. See Grounding Instructions.
10. No modification of this equipment is allowed.
11. Do not modify this equipment without authorization of the manufacturer.
12. Mattress covers have passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.

CAUTION -

1. If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.

NOTE, CAUTION AND WARNING STATEMENTS:

NOTE - Indicate some tips.

CAUTION - Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

WARNING - Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.


1. Introduction

This manual should be used for initial set up of the system and for reference purposes.

1.1 GENERAL INFORMATION

The system is a lateral rotation and true low-air-loss mattress system suitable for stage IV, bed sore treatment and prevention up to stage III.

The system has been tested and successfully approved with the following standards:

EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the IEC/EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation distances between the equipments.
- Connect the equipment into an outlet on a circuit different from that to which the other device is connected.
- Consult the manufacturer or field service technician for help.

1.2 INTENDED USE

This product is intended to help with prevention and treatment of pulmonary and skin complications related to immobility. It also provide following purposes:

- Individual acute care setting and long-term care of whom suffering from pressure ulcer.
 - Pain management as prescribed by physician.
 - To be operated by personnel who is qualified to perform general nursing procedures and has receive adequate training in knowledge of prevention and treatment of pressure ulcer.

 **NOTE:** Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

 **WARNING:** Do not rotate patients with unstable fractures, patients with unstable spinal cord injuries, or long bone traction.

 **CAUTION:** Equipment should be applied with cautions for patients with unstable internal cranial pressure or ARDS.


SYMBOLS

	Authorised representative in the European community
IPX0	ordinary use – no water ingress tested.
	Manufacturer
	Refer to instruction manual/ booklet
	"BF" symbol, indicate this product is according to the degree of protecting against electric shock for type BF equipment
	Grounding terminal
	Attention, should read the instructions!
	Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

2. PRODUCT DESCRIPTION

2.1 Pump & Mattress System

1. Quick connector
2. Pump unit
3. Turning Mattress


2.2 Pump Unit

Front


1. Quick Connector
2. Front Panel
3. Handle

Rear


1. Power Switch
2. Power Cord Socket
3. Fuses
4. Mounting Bracket
5. Pump bumper
6. Air Filter

2.3 FRONT PANEL

1. Comfort Level

Comfort level controls the air pressure output. When turning up (), the output pressure will increase and higher pressure output will support the heavier weight user. Vice versa () for decreasing air pressure. Check to see if the suitable pressure is selected by sliding one hand between the air cells and the patient to feel patient's buttocks. Users should be able to feel the minimum contact. Always leave at least 1-inch space between user's buttock areas and air cells under to prevent bottoming out.

2. Panel Lock

Should the panel remain untouched for 5 minutes, the panel lock feature will lock the panel with green LED light on to prevent accident from changing the setting during normal operation. To unlock, simply by pressing the Panel Lock button for 3 seconds.

3. Max Firm

This mode will be forced to execute every time the power is turned on. This insures the pump is able to reach its maximum operating pressure. Once the max pressure level is reached, the pump will automatically switch into the previous setting. User can also use this function as full mattress inflation during patient ingress/egress for better support.

4. On/off

Press On to turn on the unit. Press Off to turn off/standby the unit.

 **NOTE:** The power switch on the side of pump must be turned on.

5. Turning

Pressing Turning Mode and Turning Angle button sets patient rotation function and angle.

5.1 Mode

Turning mode provides 5 different settings – Pulsation, Static, Left only, Right only and both. To choose the Pulsation, press mode key several times until Left Turn, Static & Right Turn indicators light up. To choose the Static, press mode key several times until Static indicator lights up and it means no rotation. To choose Left Turn/Right Turn only, press mode key several times until Left Turn/Right Turn indicator lights up. To choose the mode of both left turn and right turn, press mode key several times until the Left Turn & Right Turn indicators light up. It results in turning from right to supine to left and then back.

5.2 Angle

Press this button sets rotation angle to half turn or full turn.

6. Dwell Time

Pressing the Dwell Time button sets the time interval that the mattress will remain in each position. Selecting either of three time intervals – 5 minutes, 30 minutes or 60 minutes – applies the same interval to all three positions - Left, Right, and Flat.

7. Alarm Mute

Press alarm mute button to temporarily suspend the CPR and Service alarm (LED light and buzzer). Should the situation not resolve within 3 mins, the alarm shall resume to notify the patient.

7.1 CPR

During the CPR situation, the caregiver can rotate the quick connector to perform CPR. Once the quick connector is switched to CPR mode, CPR alarm will sound with LED light on to notify the caregiver.

7.2 PFA (Power Failure Alarm)

During power failure situation, the Power Failure LED light will light on with buzzer. By pressing the mute button to disable both buzzer and LED.

7.3 LPA (Low Pressure Alarm)

During the low pressure situation (the pressure inside of air mattress is abnormal), the Service LED continuously lights up with buzzer (bb-bb). By pressing the mute button to disable both buzzer and LED.

7.4 Service Indicator


When Service LED continuously lights up with buzzer (b-b), the pressure inside of air mattress is abnormal or the pump has mechanical failure situation. Please check all connections are connected properly and that they are correctly.

installed as per installation instructions. Should the situation not resolve, please contact local dealer for repair immediately.

3. Installation

Unpacking the box to inspect for any damage, which may have occurred during shipment. If there are any damages, please contact your dealer immediately.

3.1 PUMP & MATTRESS INSTALLATION

1. Place the mattress on top of the bed frame. Please note for the foot end.
2. Hang the pump onto bed rail (foot-end), and adjust hangers to best upright position of the pump.

⚠ NOTE: Do not place the pump on flat surface, this might block the air intake of pump unit.

3. Connect air hose connectors from air mattress to the pump unit. Make sure the connector is in the right position as the inserted diagram below. Rotate the connector clockwise to the Lock / on position. When a "click" sound is felt or heard, the connection is completed and secured.

⚠ NOTE: The pump is only operable at "Lock/On" position.

4. Plug the power cord into electrical outlet.

NOTE:

1. Make sure the pump unit is suitable for the local power voltage.
2. The plug is also served to disconnect the device. Do not position the equipment so that it is difficult to disconnect the device.

CAUTION: The pump can only be applied to the mattress recommended by the manufacturer. Do not use it for any other purpose. (applied part: air mattress)

5. Then turn the main power switch to ON position.

Several installation tips are listed below:

After installation, the extra length of the power cord, if any, should be neatly arranged to avoid any tripping accidents. The EQUIPMENT should be firmly placed at position where users/doctors can access easily.


Grounding:

Before any connection to the output connectors is made, the unit shall be connected to a protective earth conductor via the three-core main cable; the mains plug shall be inserted only into a socket outlet provided with a protective earth contact. The protective action shall not be negated by the user of an extension cord without protective conductor.

4. Operation

⚠ NOTE: Always read the operating instruction before use.

4.1 GENERAL OPERATION


1. Switch on the main power switch found from the side of the pump and press  on the control panel to turn on the power.
2. The system will automatically go into “Max Firm” mode for a few minutes of inflation. It is highly recommended not to have patient lying on the mattress during the initial inflation. The system will be performed auto-tuning function during inflation period.

⚠ NOTE: The mattress is equipped with head pillow function, therefore the head section of the mattress and side bolster will be inflated quicker than turning air cells.

3. The Service indicator (orange LED) will light up when the mattress is not fully inflated.
4. According to the weight and height of the patient, adjust the pressure setting to the most suitable level without bottoming out.
5. Set Turning Function
 - 5.1 Set Mode – mode sequence is Left Turn, Static, Right Turn, Full Turn (both Left and Right Turn ON) and Pulsation (Static, Left and Right Turn indicators ON)
 - 5.2 Set Angle – Half Turn or Full Turn at each side.
 - 5.3 Set Dwell Time – 5, 30 or 60 minutes of pause time at each preset angle.

⚠ NOTE: Turning angle will vary depending upon the patient’s weight and lying position.

4.2 CPR

When CPR needs to be performed, quickly rotate the Quick Connector to the CPR position to release air.

4.3 PRESSURE SET UP

Users can adjust pressure of air mattress to a desired softness by adjusting the Comfort Level. Please consult your physician for a suitable setting.

⚠ NOTE: Check to see if the suitable pressure is selected by sliding one hand between the air cells and the patient to feel patient's buttocks in the static setting. Users should be able to feel the minimum contact.

4.4 TURNING FUNCTION SET UP

Users can adjust the turning mode on left and right side independently based on their personal preference and health condition. The maximum turning angle is about 40 degree (this angle may vary slightly depending on the patient's lying position and weight). The dwell time at the flat level and turn angle as 5, 30 or 60 minutes can be adjusted. Please consult your physician for a suitable setting.

For example, if the turning mode is set at left turn only with 5 min dwell time. The mattress will turn left to the preset turning-angle position and stay there for 5 min, rotate back to flat position and stay there for 5 min. This cycle will keep repeating until the setting is changed. One more example, if the turning mode is set at Full Turn (both left and right turn ON) with 30 min dwell time. The mattress will turn left first to the preset turning-angle position for 30 min, rotate back to flat position and stay there for 30 min, turn right to the same preset turning-angle position for 30 min, and then turn back to the flat position and stay there for another 30 min.

4.5 AUDIABLE AND VISIBLE ALARM

When the following situation occurred, the LED lights will light up along with the buzzer.

4.5.1 CPR – When CPR or quick storage is performed, or the Quick Connector is on Disconnect position under power on situation, the “CPR” indicator (orange LED)

will light up along with buzzer.

- 4.5.2 Power Failure – When electrical shortage occurred or power cord is disconnected without turning off the pump first, the “power failure” indicator (orange LED) will light up along with buzzer. Check to ensure power cord is connected properly.
- 4.5.3 Low Pressure – During the low pressure situation (the pressure inside of air mattress is abnormal), the Service LED continuously lights up with buzzer (bb-bb). By pressing the mute button to disable both buzzer and LED
- 4.5.4 Service – When pump has mechanical failure situation, the “Service” indicator (orange LED) will light up along with buzzer (b-b). Check to ensure all connections are connected properly and that they are correctly installed as per installation instructions. The Service indicator will continuously light up until the abnormal condition is resolved. Should the situation not resolve, please contact local dealer for repair immediately.

 **NOTE:** If the low-pressure situation is not resolved within 1 minute, the system will shutdown for safety consideration. Please check for any leakage (tubes or connecting hoses). If necessary, replace any damaged tubes or hoses. Or contact local qualified dealer for repair.

4.6 ALARM MUTE

When these situations mentioned above occurred, both the LED light and buzzer will sound to warn the patient. By pressing the button, it will temporary mute the buzzer (CPR & Low Pressure) so the caregiver may check for possible causes. If the CPR & low pressure situation not resolved within 3 minutes, the alarm will resume to notice the user. When in Power Failure situation, pressing alarm mute button, it will cease the buzzer and press the ON/OFF key to turn off the “power failure” indicator.

5. Cleaning

By wiping the pump unit with a damp cloth pre-soaked with a mild detergent, and keep it away from dust. If other detergent is used, choose one that will have no chemical effects on the surface of the plastic case of the pump unit.

 **CAUTION:** Do not immerse or soak pump unit.

By wiping the mattress unit with a damp cloth pre-soaked with warm water (not to exceed 65 °C) containing a mild detergent, and keep it away from dust. The cover may also be cleaned using sodium hypochlorite diluted in water. All parts should be air dried thoroughly before use.

 **CAUTION:** Do not use phenolic-based product for cleaning.

 **CAUTION:** After cleaning, dry the mattress without direct exposure of sunlight.

The carrying bag (if available) should be turned inside out and completely wiped down using the disinfectant solutions. Allow it to air-dry thoroughly. Once the inside is dry, turn it back and wipe down the outside of the bag with disinfectant solutions.

6. Storage

1. To quickly vacuum air out from mattress for storage, rotate the Quick Connector to the CPR position and press the "Max Firm" button on the control panel. Air inside of mattress will be vacuumed out within 5 minutes and then pump will be stopped automatically.
2. Rotate Quick Connector to the "Disconnect" position and pull up the Quick Connector.
3. Lay the mattress out flat and upsides down.
4. Roll from the head end towards the foot end.
5. Foot-end strap can then be stretched around the rolled mattress to prevent unrolling.
6. The power cord could be wrapped around the pump bumper on the back of pump.

 **CAUTION:** When the pump is vacuuming air, air temperature is hot. Do not block the air outlet.


 **NOTE:** Do not fold, crease or stack the mattresses.

7. Maintenance

7.1 GENERAL


1. Check main power cord and plug, if there are abrasions or excessive wears.
2. Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are stubbed together correctly.
3. Check the air hoses for any kink or break. For replacement, please contact your local dealers.

7.2 FUSE REPLACEMENT

1. Disconnect the plug from mains power when a blown fuse is suspected.
2. Remove the cover of the fuse holder by means of a small screwdriver.
3. Insert a new fuse of the correct rating in, and replace the cover of the fuse holder back. The fuse should be rated as T3.15AL /250V (for 230V system) and T5A /250V (for 120V system) type .

7.3 AIR FILTER REPLACEMENT

1. Replace the air filter located at the side of pump.
2. The filter is reusable and can be washed gently with a mild detergent and water. Dry the filter before use.
3. Check and replace air filter regularly if environment is dirty.

7.4 RECHARGEABLE BATTERY

1. To check if the rechargeable battery has been drained out, unplug the power cord and see if the Power Failure indicator will light up along with buzzer last for a few minutes.
2. If the Power Failure Alarm is unable to work or the battery might need to be replaced (approximate life expectancy 6 months), user can contact your dealer or notify the technician for replacement.

8. TROUBLESHOOTING

Q1.Power is not ON

- Check if the plug is connected to mains.
- Check for a blown fuse.

Q2.Alarm is on(Audible & visual)

- Check if the Quick Connector is rotated to the Lock/on position.
- Check if the power is suddenly shut down.
- Check if the connection between air tube connector to pump unit is tightly secured.
- Check if all tubing connections along mattress are secured.

Q3.Patient is bottoming out

- Pressure setting might be inadequate for the patient, adjust comfort range 1 to 2 levels higher and wait for a few minutes for best comfort.

Q4.Mattress seems not turning

Wait for a complete turning cycle. It might be at the stage of flat level or turn level dwell interval.

If the above information does not solve your problems, please contact your local agent directly. They might require a technician to take care the problem.

9. EXPECTED SERVICE LIFE

The products are intended to offer safe and reliable operation when use or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indications on products. Otherwise, service and inspection of the devices generally should not be required.

10. Technical Description

PUMP		SPECIFICATION
Power Supply (Note: See rating label on the product)		AC 100-120V 50/60 Hz, 5.0A (for 120V system) AC 220-240V 50/60Hz, 3A (for 230V system)
Fuse Rating		T5A /250V (for 120V system) T3.15AL /250V (for 230V system)
Dimension (L x W x H)		15.8 x 38 x 26.8 (cm) or 6.2" x 15" x 10.5"
Weight		6 kg or 13.2 lb
Classification		Class I, Type BF, IPX0 Applied Part: Air Mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)
ENVIRONMENTAL INFORMATION		
Environment	Operating Altitude	2000 m (Maximum)
	Temperature	Operation: 10° C to 33° C (50° F to 91° F) Storage: -15° C to 50° C (5° F to 122° F) Shipping: -15° C to 70° C (5° F to 158° F)
	Humidity	Operation: 10% to 90% non-condensing Storage: 10% to 90% non-condensing Shipping: 10 % to 90% non-condensing
MATTRESS		SPECIFICATION
Model		8" Mattress with 10" side bolster and 3 drainage sites 2" foam pocket is provided
Dimension (L x W x H)		200 x 90 x 25.4 cm or 78.7" x 35.4" x 10"
Weight		10.1 kg or 22.2 lb
Cover		Stretch
Cell Material		PU Coated Nylon
Base Material		PU Coated Nylon

NOTE:

- 1.Consult the distributor or EU representative for further technical documents.
- 2.The specification is also suitable for other areas operating with same power supply.
- 3.Mattress dimension and weight is measured without foam cushion
- 4.The manufacturer reserves the right to modify the specification without notice.

Appendix A: EMC Information

GUIDANCE AND MANUFACTURER’S DECLARATION- ELECTROMAGNETIC EMISSIONS:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	
<p>⚠ Warning:</p> <ol style="list-style-type: none"> 1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used. 2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. 3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. 		

GUIDANCE AND MANUFACTURER'S DECLARATION- ELECTROMAGNETIC IMMUNITY:

This device is intended for use in the electromagnetic environment specified below.

The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	±8kV contact ±15kV air		±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply line ±1kV for input/output line		±2kV for power supply line ±1kV for input/output line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25/30 cycles <5 % UT (>95 % dip in UT) for 250/300 cycles		120V/240V	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz)	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of

<p>magnetic field IEC61000-4-8</p>				<p>a typical location in a typical commercial or hospital environment.</p>
<p>Conducted RF IEC 61000-4-6</p>	<p>3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz</p>	<p>3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz</p>	<p>6Vrms</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.</p>
<p>Radiated RF EM Fields IEC61000-4-3</p>	<p>3 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz</p>	<p>10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz</p>	<p>3V/m</p>	<p>Recommended separation distance $d = 2.0 \sqrt{P}$ 150kHz to 80MHz $d = 2.0 \sqrt{P}$ 80MHz to 800MHz $d = 4.0 \sqrt{P}$ 800 MHz to 2.7GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency ranged.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1: U_T is the a.c. mains voltage prior to the application of the test level
NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 2.0\sqrt{P}$	80 MHz to 800 MHz $d = 2.0\sqrt{P}$	800 MHz to 2.7 GHz $d = 4.0\sqrt{P}$
0.01	0.2	0.2	0.4
0.1	0.63	0.63	1.26
1	2	1	4
10	6.3	6.3	12.6
100	20	20	40

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Print-2018/All rights reserved
Apex Medical Corp.