
WellNeb User Manual

Be well, Live well
Wellell

TABLE OF CONTENTS

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PLEASE READ THE MANUAL BEFORE USE

IMPORTANT SAFEGUARDS

Please make this precaution list available to all users and providers

**READ ALL INSTRUCTIONS BEFORE USE NOTE, CAUTION, AND
WARNING STATEMENTS**

DANGER: To reduce the risk of electrocution

- (1) DO NOT place or store product where it can fall into or be pulled into a tub or sink.
- (2) DO NOT place or immerse the unit into water or other liquids.
- (3) Do not use while bathing.
- (4) Do not reach for a product that has fallen into water. Unplug immediately.

WARNINGS: To reduce the risk of burns, electrocution, fire, or injury to the user or operator.

- (1) Use this product only for its intended purpose as described in this manual.
- (2) The user or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or patient is established.
- (3) Never perform servicing or maintenance to the device while it is in use.
- (4) Never disassemble this device on your own. Repairs and internal servicing should only be performed by an authorized technician.
- (5) THIS DEVICE IS NOT INTENDED FOR LIFE SUPPORT.
- (6) Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with nitrous oxide or oxygen.
- (7) Do not leave unit unattended while plugged in.
- (8) If any abnormality occurs, discontinue use immediately until the unit has been examined and repaired.
- (9) Always unplug the product immediately after use.
- (10) Never block the air openings of the main unit or place it where the air openings may be obstructed.
- (11) Do not operate the unit with an empty medicine cup.

- (12) Adult supervision is highly recommended when used by children.
- (13) Do not direct medication mist output toward eyes.
- (14) Do not use while operating a vehicle.
- (15) Do not overfill the medicine cup.
- (16) Always comply with the local regulations for disposing of such equipment.
- (17) Store the product in a dry place. Keep out of direct sunlight. Use of or store the product only in the following conditions:
 - Operation: 5°C to 40°C (41°F to 104°F); 15%RH to 90%RH, non-condensing
 - Storage: -25°C to 70°C (-13°F to 158°F); 0% RH to 90% RH, non-condensing

CAUTION: Reducing the risk of damage to the device.

- (1) To avoid electrical shock: Keep the unit away from water.
- (2) Do not immerse the power cord or the unit in liquid. Do not use while bathing.
- (3) Do not reach for a unit that has fallen into water – immediately unplug the unit.
- (4) Never operate the unit if it has any damaged parts (including power cord) or if it has been dropped or submersed in water. Promptly send it to a service center for examination and repair.
- (5) The unit should not be used where flammable gas, oxygen or aerosol spray products are being used.
- (6) Disconnect the unit from the electrical outlet before cleaning, filling, and after each use.
- (7) Do not use the other attachments unless recommended by manufacturer.

NOTE: Indicates information that you should pay special attention to.

-  **CAUTION: Indicates correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.**
-  **WARNING: Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.**
-  **Caution: Federal law (USA) restricts this device to sale by, or on the order of, a licensed physician. Carefully read instructions before use.**

SYMBOLS

	CE marking		Medical Device
	Authorized representative in the European community		Type BF applied parts: mouthpiece, masks
IP21	Protection against harmful ingress of water: The product is not protected against any physical contact or objects.		Class II Equipment
	Manufacturer		Date of manufacture
	Refer to instruction manual or booklet		Serial number
	This way up		Fragile, handle with care
	General symbol for recovery recyclable		Use no hand hooks
	Humidity limitation		keep away from sunlight
	Keep dry		Caution
	To indicate the entity importing the medical device into the locale		
	<p>Attention</p> <p>Observe proper Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.</p>		

1. Introduction

The device should be used under the supervision and direction of a licensed physician and/or a respiratory therapist.

Thank you for purchasing the WellNeb. With proper care and use, it will provide you with many years of reliable treatment.

This product is developed for the successful treatment of asthma, allergies and other respiratory disorders. It creates a stream of air that travels through the clear tube to the nebulizer. When air enters the nebulizer, it will convert the prescribed medication into an aerosol mist for easy inhalation.

We encourage you to read this guidebook thoroughly to learn about the features of this product. Any use of this product other than its intended use should be avoided.

1.1. Intended Use

Your Nebulizer is intended for use in the treatment of respiratory ailments in which aerosolized medication is required. Please consult your physician to determine if your prescription medication is approved for use with nebulizers.

Single Patient Use Only.

2. Product Description

2.1. Content of Package

Unpack to check for any damage which may have occurred during shipment. If there are damages, please contact your dealer immediately.

This system comes with the following items.


- (1) Pump compressor unit x1
- (2) Air tube x1 (7ft)
- (3) Nebulizer
- (4) Angled mouthpiece
- (5) Air filter x5
- (6) User Manual x1

 **Polarization Instructions:** This device has a polarized plug (one blade is wider than the other). To reduce the risk of electric shock, this plug is intended to fit into a polarized outlet only one way. If the plug does not fit fully into the outlet, reverse the plug. If it still does not fit, contact a qualified electrician.

 **Do not attempt to modify the plug in any way.**

 **Replacement parts are available through your local retailer. Optional accessories include adult and child masks.**

2.2. Product Overview

- (1) Filter and Cover
- (2) Tube Connector
- (3) Power Switch
- (4) Nebulizer


- (5) Air Tube
- (6) Mouthpiece
- (7) Mask (Adult/child)

3. Before Use

 **Note: Prior to initial operation, the nebulizer should be cleaned thoroughly. Refer to the “Cleaning” section in the manual.**

- (1). Place your compressor nebulizer on a flat and stable surface. Be sure that you can easily reach the controls when you are seated.
- (2). Get accessories ready.
- (3). Gently twist the top part of the nebulizer counterclockwise to disassemble the nebulizer.
- (4). Fill the bottom section of the nebulizer with the medication prescribed by your physician. Be sure that the purple cone is put inside the bottom section as shown.
- (5). Gently twist the top part clockwise to reassemble the nebulizer. Be sure the two sections fit well.
- (6). Attach one end of the air tube to the base of the nebulizer.
- (7). Attach the other end of the air tube to the air tube connector located on the back of the compressor.
- (8). Attach mouthpiece to the top section of the nebulizer. If using a mask, attach the nebulizer bottle to the opening in the bottom of the mask.
- (9). Plug the power cord into an appropriate electrical outlet. Make sure at this stage, the POWER switch is at the “off” status.
- (10). Press POWER switch to begin your prescribed treatment. Important: The compressor motor has a thermal protector that will shut off the unit before the unit is overheated. When the thermal protector shuts the unit off, please:
 - Switch off the unit.
 - Unplug the unit from the electrical outlet.
 - Wait 30 minutes for the motor to cool down before another treatment. Make sure the air openings are not obstructed.
- (11). When treatment is finished, shut off the unit and unplug it from the electrical outlet.


 **Place purple cone inside and fill medication.**

4. Operation

 **Note: The nebulizer is operable up to a 45° angle.**

- (1). Attach one end of the air tube to the air tube connector.
- (2). Attach the opposite end of the air tube to the stem at the base of nebulizer.
- (3). Attach the angled mouthpiece or optional mask to the top of the nebulizer. If using an optional tee adapter set:
 - Attach the tee adapter to the top of the nebulizer.
 - Connect one end of the tee adapter to the straight mouthpiece and the other to the extension tube.



 **Note: The tee adapter is not needed if using a mask.**

- (4). Check that the POWER switch is in the “off” position. Plug power cord into outlet.
- (5). Press the POWER switch to the “on” position.
- (6). Start your treatment according to your physician’s or therapist’s instruction.

 **Note: The small groove on the handle of the compressor can be used to support the nebulizer before, during or after treatment.**

- (7). When the treatment is complete, press the POWER switch to the “off” position and unplug the power cord from the outlet.


5. Cleaning

It is recommended that the nebulizer and mouthpiece are thoroughly cleaned with hot water after each use and cleaned with a mild detergent after the last treatment of the day.

If your physician or respiratory therapist specifies a different cleaning procedure, follow their instructions.

- ⚠ **Check all unit components (air filter, nebulizer, mouthpiece, mask) before each use.**
- ⚠ **Do not immerse the unit in water. It may damage the unit.**
- ⚠ **Disconnect the unit from the electrical outlet before cleaning.**
- ⚠ **Clean all necessary parts after each use as instructed in this guidebook.**

5.1. Rinsing (after each use)

- (1). Make sure that the power switch is turned off and that the unit is unplugged.
- (2). Disconnect the air tube, nebulizer and mouthpiece.
- (3). Gently twist the nebulizer to open it.
- (4). Rinse the nebulizer and mouthpiece with hot water.
- (5). Dry with clean, soft towel or let it air dry.
- (6). Reassemble the nebulizer when completely dry and put these parts in a dry, sealed container.

5.2. Cleaning (the last treatment of the day)

Please follow these steps to disinfect your nebulizer unless otherwise specified by your physician. It is suggested that the unit is disinfected after the last treatment of the day.


- (1). Follow steps 1-3 under "Rinsing."
- (2). Wash the nebulizer parts and mouthpiece in warm water and a mild detergent. Then, wash them in hot tap water.
- (3). Use one part white vinegar with three parts distilled water. Make sure the mixed solution is enough to submerge the nebulizer and mouthpiece.
- (4). Submerge these parts in the vinegar and water solution for 30 minutes.
- (5). Complete rinsing steps 4-6.

- ⚠ **Cleaning the Compressor: Wipe daily with a damp cloth.**
- ⚠ **Note: Do not use any powdered cleaners or soap pads, which may damage the finish.**

6. Assembly

Follow the cleaning instructions in this manual under “Cleaning” prior to use for the first time, or after unit has been stored for an extended period of time.

- (1). The compressor should be on a flat stable surface and within reach.
- (2). Gently twist and pull straight up on the lid of the nebulizer to separate the medicine cup and cover.
- (3). Make sure the purple cone is properly installed inside the medicine cup. The stem inside the medicine cup inserts into the tube of the purple cone.
- (4). Add the prescribed amount of medication to the medicine cup.
- (5). Reassemble the nebulizer by carefully twisting the medicine cup and cover together. Both parts should fit securely.

7. Maintenance

7.1. Storage

- Keep away from sunlight. Keep dry.
- Keep the unit out of reach of small children.
- Always keep the unit unplugged while not in use.

7.2. Filter Replacement

It should be replaced every 30 days of use, or sooner if it appears dirty.

- (1). Open the blue filter cover.
- (2). Discard the used filter.
- (3). Replace the used filter with a new, clean filter.
- (4). Close the filter cover.

WARNING: The filter is NOT washable.

WARNING: Use only the manufacturer authorized filter as a replacement.

⚠ CAUTION: Dirty air filter may cause high operating temperatures that might impair device performance.

8. Specifications

Item	Specification
Power	230V, 50Hz
Power Consumption	180VA
Medication Capacity	6mL
Particle Size	0.5 µm to 10 µm
MMAD	5µm
Sound Level	Less than 60 dB at 1 meter
Average Nebulization Range	≥0.2 mL/min
Compressor Pressure Range	≥25 psi
Operating Pressure Range	8 to 13 psi
Liter Flow Range	5.0 Lpm ~ 8.0 Lpm
Operating Temperature Range	5°C to 40°C (41°F to 104°F)
Operating Humidity Range	15 RH - 90% RH
Storage Temperature Range	-25°C to 70°C (-13° to 158°F)
Storage Humidity Range	0% RH - 90% RH
Dimension (L x W x H)	238 mm x 178 mm x 99 mm (9.37" x 7.01" x 3.90")
Weight	1.36 kgs
Standard Accessories	Nebulizer, air tube, mouthpiece, filters Optional Accessories Adult mask, child mask

9. Disposal

Observe proper Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

Always comply with the local regulations for disposing of such equipment.

10. Servicing

The products are intended to offer safe and reliable operation when use or installed according to the instructions provided by Wellell. Wellell recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indication on products. Otherwise, service and inspection of the devices generally should not be required.

11. Limited Warranty

In the case of a complaint, **the consumer should contact the commercial establishment where the product was purchased or the authorised supplier.**

In the event that it is impossible to contact the establishment where the product was purchased, for any queries or complaints you can send an e-mail to incidencias@wellell.com, call **94.470.64.08** or contact our sales representatives.

The guarantee of our products starts from the day WELLELL delivers the product to the customer, generally from the date of the delivery note plus 1 or 2 days for transport.

The product will be collected or delivered, if necessary, at the customer's place where the corresponding invoice has been issued. The product must be correctly packaged and in a reasonable hygienic condition.

For the approval of the warranty as a result of a complaint or claim for any defect or incident in the product, **the following information must be submitted:**

- The name or company name of the purchaser.
- The serial number of the product.
- Description of the problem or defect of the product.
- Photos of the damaged equipment in the event of breakage, damage, knocks, etc.

Failure to provide this information will exempt the warranty.

Defects or damage due to incorrect storage, use or handling of the product or its materials or wear and tear caused by normal use of the product are not included in this warranty.

Wellell disclaims any liability for personal injury or property damage resulting from improper or dangerous use of this chair, lack of maintenance or failure to follow the instructions in the user manual.

The purchaser assumes and accepts all conditions of sale and warranty set forth herein in the purchase of each product purchased from WELLELL.

Information of Electromagnetic compatibility (EMC)

Guidance and manufacturer's declaration –electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations	<5 % UT (>95 % dip in UT) for 0.5 cycle	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the

on power supply input lines IEC 61000-4-11	40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s		user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from a continuous power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2\sqrt{P}$ 80MHz to 800MHz $d = 2.3\sqrt{P}$ 800MHz to 2.5GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80MHz and 800MHz, the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d = 1.2 \sqrt{P}$	80MHz to 800MHz $d = 1.2\sqrt{P}$	800MHz to 2.5GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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All rights reserved.
Biolife Ltd.
Vega Technologies Inc.