Domus 2+

User Manual

Be well, Live well

Wellell

CONTENTS

ALTERNATING PRESSURE REDISTRIBUTION SYSTEM	
USER MANUAL	1
SISTEMA DE REDISTRIBUCIÓN DE PRESIÓN ALTERNA	
MANUAL DE INSTRUCCIONES	18
SISTEMA DE REDISTRIBUIÇÃO DE PRESSÃO ALTERNADA	
INSTRUÇÕES DE UTILIZAÇÃO	35
SYSTÈME DE REDISTRIBUTION À PRESSION ALTERNÉE	
INSTRUCTIONS D'UTILISATION	52
ALTERNIERENDES DRUCKUMVERTEILUNGSSYSTEM	
GEBRAUCHSANLEITUNG	69
ALTERNEREND DRUKSYSTEEM	
GEbruikersHANDLEIDING	87
SISTEMA DI RIDISTRIBUZIONE A PRESSIONE ALTERNATA	
MANUALE DI ISTRUZIONI	105

MODEL NO.: 9P-046520

PLEASE READ THE MANUAL BEFORE USE

IMPORTANT SAFEGUARDS

SAVE THESE INSTRUCTIONS

READ ALL INSTRUCTIONS BEFORE USING

DANGER - To reduce the risk of electrocution:

- (1) Always unplug this product immediately after using.
- (2) Do not use while bathing.
- (3) Do not place or store this product where it can fall or be pulled into a tub or sink.
- (4) Do not place in or drop into water or other liquid.
- (5) Do not reach for a product that has fallen into water. Unplug immediately.

WARNING - To reduce the risk of burns, electrocution, fire, or injury to persons:

- (1) Evaluate patients for entrapment risk according to protocol and monitor patients appropriately.
- (2) This system is not for use with patients who have a spinal cord injury.
- (3) Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may result from a child swallowing a small part detached from the device.
- (4) Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer
- (5) Never operate this product if it has a damaged cord or plug, if it is not working properly, or if it has been dropped, damaged, or dropped into water. Return the product to a service center for examination and repair.
- (6) Keep the cord away from heated surfaces.
- (7) Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where their openings may be blocked. Keep the air openings free of lint, hairs, and other similar particles.
- (8) Never drop or insert any object into any opening or hose.
- (9) Do not modify this equipment without authorization of the manufacturer.
- (10) The covers have passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.
- (11) Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.

CAUTION –

If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.

NOTE: Indicate some tips.

CAUTION: Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property

WARNING: Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

SYMBOLS

Authorized representative in the European community

Catalog, reorder or reference number

Manufacturer

Complies with standards protecting against electric shock for type BF equipment.

Attention, you should read the accompanying information carefully!

Attention – Observe proper Disposal of Electrical & Electronic Equipment (WEEE):

This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

Consult operating instructions for use

IP2X

Protected against solid foreign objects of 12,5 mm and greater; No protection against vertically falling water drops

Class II

Temperature limitation/temperature range

Do Not Iron

Tumble Dry, Normal, Low Heat

Do Not Tumble Dry

Do Not Bleach

Do Not Dry Clean

Machine wash, regular / normal, 60 degrees C (140 degrees F)

KEEP THIS MANUAL FOR REFERENCE

1. Introduction

This manual should be used for initial set up of the system and for reference purposes.

1.1. General Information

The system is a high quality and affordable mattress system suitable for treatment and prevention of pressure ulcers.

The system has been tested and successfully approved to the following standards:

IEC/EN 60601-1
IEC/EN 60601-1-2
IEC/EN 61000-3-2 Class A
IEC/EN 61000-3-3
CISPR 11 Group 1, Class B

EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

1.2. Intended Use

This product is intended to:

- to help and reduce the incidence of pressure ulcers while optimizing patient comfort.
- for long term home care of patients suffering from pressure ulcers.
- for pain management as prescribed by a physician.

The product can only be operated by personnels who are qualified to perform general nursing procedures and has received adequate training in knowledge of prevention and treatment of pressure ulcer.


NOTE: Equipment not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.

2. Product Description


Unpack the box to check for any damage which may have occurred during shipment. If there are damages, please contact your dealer immediately.

2.1. Pump and Mattress System

- (1) CPR*
- (2) Pump unit
- (3) Mattress unit
- (4) Quick Connector


2.2. Pump Unit

Front


- (1) Quick Connector Slot
- (2) Display Panel
- (3) Quick Connector with Cap

Rear

- (4) Power Cord
- (5) Clip Hanger
- (6) Fuses
- (7) Filter

2.3. Front panel

(1) pressure Adjustment Knob

The pressure adjustment knob controls the air pressure output. When turning clockwise, the output pressure will be increased. On the other hand, when turning counter-clockwise, the output air pressure will be decreased.

⚠ NOTE: You can also check if the pressure is suitable for the patient by sliding one hand beneath the air cells at the level of the patient's buttocks. The air cells will alternately inflate and deflate. You should feel slight contact with the buttocks when the air cells beneath the buttocks deflate.

(2) Low Pressure Indicator

When the mattress is initializing, the Low Pressure LED will illuminate until the appropriate pressure is reached (according to the Pressure Adjustment Knob); this is normal. Otherwise, the Low Pressure LED is a warning, indicating that the pressure in mattress is unusually low. Check to ensure all connections are secured and correctly installed as per the installation instructions.

(3) Normal Pressure Indicator

When the green LED illuminates, the pressure inside of air mattress has reached the desired pressure setting (preset on the Pressure Adjustment Knob).


(4) Main Power Switch

To turn on/off the pump unit with illuminated green light.

3. Installation

3.1. Pump & Mattress Installation

Several installation tips are listed below:


- (1) Place the mattress on top of the bed frame. Fasten the Nylon Straps to secure the mattress to bed frame if available.


Warning: The overlay mattress must be applied on the underlying mattress

- (2) Clip the plastic hanger on the back of pump in a bottom-to-top direction and hang the pump onto bed rail (foot-end) or place the pump on a flat surface.
- (3) Connect quick connector from air mattress to the pump unit. When a "click" is felt or heard, the connection is completed and secured.

NOTE: Check and ensure the air hoses are not kinked or tucked under mattress.



Warning: After installation, the extra length of the power cord, if any, should be neatly arranged to avoid any tripping accidents. The EQUIPMENT should be firmly placed at position where users/doctors can easily access.

- (4) Plug the power cord into electrical outlet.

 **NOTE:**

1. Make sure the pump unit is suitable for the local power voltage.

2. The plug is also served as the disconnect device. Do not position the equipment so that it is difficult to operate the disconnecting device.


 **CAUTION:** The pump can only be applied to the mattress recommended by the manufacturer. Do not use it for any other purpose. (applied part: air mattress)

- (5) Then turn the main power switch to ON position which located on the display panel.

4. Operation

 **NOTE: Always read the operating instruction before use.**

4.1. General operation


- (1) Switch on the main power switch on the display panel.

- (2) The pump unit is beginning to deliver air into mattress. Turn the comfort knob clock-wise to increase firmness. (i.e. pressure)


- (3) The low pressure indicator (yellow LED) will light up when the mattress is not yet fully inflated.

- (4) When the appropriate pressure is reached, the low pressure indicator (yellow LED) will go off and normal pressure indicator (green LED) will light up automatically.

 **NOTE: Each time when mattress is firstly inflated, it is recommended that the pressure knob be set to "Firm" for the quickest inflation. User can then adjust air mattress to the desired softness afterward.**

4.2. Transport

For transport purpose, interconnect Male to Female part of the air hose connector. When a "click" is felt or heard, the connection is completed and secured; then air from mattress is sealed off.

5. Cleaning

It is important to follow the cleaning procedures before using the equipment on human bodies; otherwise, patients and/or doctors may have the possibility of getting infection.

Wipe the pump unit with a damp cloth and a mild detergent, and keep it away from dust. If other detergent is used, choose one that will have no chemical effects on the surface of the plastic case of the pump unit.

 **CAUTION: Do not immerse or soak pump unit.**

Wipe down the mattress with warm water containing a mild detergent. The cover may also be cleaned by using sodium hypochlorite diluted in water. All parts should be air dried thoroughly before use.

Cover Material: Nylon/PU	    
Cover Material: Silver+ Nylon/PU	

 **CAUTION: Do not use phenolic based product for cleaning.**

 **CAUTION: Dry the mattress on SUNLESS capitalized area after cleaning.**

6. Storage

- (1) To store the mattress, lay the mattress out flat and upsides down.
- (2) Roll from the head end towards the foot end with CPR valve open.
- (3) Foot-end strap can then be stretched around the rolled mattress to prevent unrolling.

 **NOTE: Do not fold, crease or stack the mattresses.**


Avoid direct sunshine.

7. Maintenance

7.1. General


- (1) Check main power cord and plug if there is abrasions or excessive wear.
- (2) Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are stubbed together correctly.
- (3) Check airflow from the quick connector. The airflow should alternate between each connector every half-cycle time if it's in alternating mode.
- (4) Check the air hoses if there is kink or breaks. For replacement, please contact our local dealers.

7.2. Fuse Replacement

- (1) Disconnect the plug from mains power when a fuse is to be replaced.
- (2) Remove the cover of the fuse holder by means of a screwdriver. (If necessary, un-clip the hanger from the pump unit)
- (3) Insert a new fuse of the correct rating in, and replace the cover of the fuse holder back. The fuse should be T1AL/250V type and VDE approved.

7.3. Air Filter Replacement

- (1) Open the air filter plate located at the back of pump.
- (2) Suggest to check and replace air filter regularly
- (3) if environment is dirty.

8. Troubleshooting

If your questions can't be answered with above information, please contact your local agent directly. They might require a technician to take care of the problem.

Problem	Solution
Power is not ON	<ul style="list-style-type: none"> • Check if the plug is connected to mains. • Check if there is any blown fuse.
Low pressure indicator is ON	<ul style="list-style-type: none"> • Check if the connection between quick connector to pump unit is tightly secured. • Check if all tubing connections along mattress are secured. • Check if the CPR valve is sealed. • Check if there are any air leakage on air cells.
Patient is bottoming out	<ul style="list-style-type: none"> • Pressure setting might be inadequate for the patient, adjust pressure to higher level and wait for a few minutes.
Mattress form is loose	<ul style="list-style-type: none"> • Check if all the snap buttons or straps of mattress are all securely fastened. • Check if the mattress is fixed to the bed frame by straps.
No air produced from some air outlets of the air tube connector	<ul style="list-style-type: none"> • This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time.

9. Expected service life

The products are intended to offer safe and reliable operation when use or installed according to the instructions provided by Wellell. Wellell recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indication on products. Otherwise, service and inspection of the devices generally should not be required.

10. Technical Description

Item		Specification
Power Supply (Note: Refer to the rating label on the product)		AC230V 50/60Hz,0.05A (for 230V system)
Fuse Rating		T1AL, 250V
Cycle time	PM100B Model: 9P-046520	8 min for 60 Hz / 9.6 min for 50 Hz
Dimension(LxWxH)		25 x 13.5 x 9.5 (cm) or 9.8" x 5.3" x 3.7"
Weight		1.8 Kg or 4.0 lb
Environment	Atmospheric Pressure	700 hPa to 1013.25 hPa
	Temperature	Operation:10 °C to 40 °C (50 °F to 104 °F) Storage: -15 °C to 50 °C (5 °F to 122 °F) Shipping: -15 °C to 70 °C (5 °F to 158 °F)
	Humidity	Operation:10% to 90% non-condensing Storage:10% to 90% non-condensing Shipping:10% to 90% non-condensing
Classification		Class II, Type BF, IP2X (for 230V system) Applied Part: Air Mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)
Mattress		Speification
Model		Domus 2+
Dimension (LxWxH)		200x80x10.2 (cm) or 78.7" x 31.5" x 4"
Weight		3.8 kg or 8.4 lb

NOTE:

- (1) **The specification is also suitable for other areas operating with same power supply.**
- (2) **Consult the distributor or EU representative for other technical documents.**
- (3) **The manufacturer reserves the right to modify the specifications without notice.**

Appendix A: EMC Information

Guidance and Manufacturer's Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	
 Warning: <ol style="list-style-type: none"> 1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used. 2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. 3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. 		

Guidance and Manufacturer’s Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	±8kV contact ±15kV air		±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply line ±1kV for input/output line		±2kV for power supply line ±1kV for input/output line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	± 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		230V	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	6Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.

<p>Radiated RF EM Fields IEC61000-4-3</p>	<p>3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz</p> <p>385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation</p>	<p>10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz</p> <p>385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation</p>	<p>10V/m</p>	<p>Recommended separation distance</p> $d = \sqrt{P}$ <p>150kHz to 80MHz</p> $d = 0,6\sqrt{P}$ <p>80MHz to 800MHz</p> $d = 1,2\sqrt{P}$ <p>800 MHz to 2.7GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1: U_T is the a.c. mains voltage prior to the application of the test level</p> <p>NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people</p>				
<p>a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.</p> <p>b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.</p>				

Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \sqrt{P}$	80 MHz to 800 MHz $d = 0,6\sqrt{P}$	800 MHz to 2.7 GHz $d = 1,2\sqrt{P}$
0.01	0.1	0.06	0.12
0.1	0.31	0.19	0.38
1	1	0.6	1.2
10	3.1	1.9	3.8
100	10	6	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Print-2024-01-22/All rights reserved
Wellel Inc.